

Instructions for: **SOIL/HW TCLP Carbamates (Cat # PT-TCLPCARB-SOIL)**

SCOPE AND APPLICATION

Wibby Environmental's SOIL/HW TCLP Carbamates Proficiency Testing Standard is designed to be extracted using USEPA Method 1311 followed by analysis for Carbamates. If you have any questions about the use of this standard, please contact Wibby Environmental Customer Service at 1-866-WibbyPT (866-942-2978).

SAMPLE PREPARATION AND ANALYSIS

1. **Preliminary evaluations (EPA Method 1311, Section 7.1).** You do not need to complete the preliminary evaluations of the Wibby Environmental TCLP Carbamates standard. The results of these tests are:
 - a. Percent solids = 100%
 - b. Particle size reduction = not needed
 - c. Required extraction fluid = fluid #1
2. For best results, the PT standard spiking concentrate should be stored at 4°C and warmed to room temperature (near 20°C) prior to use.
3. The SOIL/HW TCLP Carbamates standard includes two separate containers, a jar containing Blank Soil Matrix and a 2mL flame sealed ampule containing the TCLP Carbamates spiking concentrate.
4. Open the Blank Soil Matrix jar and transfer 100 g of the matrix to your TCLP extraction vessel.
5. Carefully break the top of the spiking concentrate ampule and add 1.0 mL (milliliters) of the concentrate using a gastight syringe onto the Blank Soil Matrix in the TCLP Extraction vessel.
6. As specified in Method 1311 add 2000 mL (20X the sample size of 100 g) of TCLP extraction fluid #1 into the extraction vessel.
7. The standard is now ready for extraction by Method 1311 and analysis per the selected method(s).

NOTE: If your normal procedure for TCLP Carbamates uses a sample size that differs from 100 g adjust the amount of Blank Soil Matrix, TCLP Carbamates spiking concentrate and TCLP Extraction Fluid as appropriate. For example, if your method calls for the extraction of a 50 g sample, use 50 g of Blank Soil Matrix, 500 µL (microliters) of TCLP Carbamates spiking concentrate and 1000 mL of TCLP Extraction Fluid #1.

REPORTING RESULTS

1. Report all results on line at www.wibby.com. Click on the "Online data entry link" and follow the instructions to enter your results. As an alternative, you may enter your results using the Data Reporting Sheets enclosed with your standards.
2. Report your results to three significant figures in ug/L of TCLP extraction fluid.
3. FAX your results to Wibby Environmental at 303-940-0043. You may also mail your results to Wibby Environmental, 6390 Joyce Drive, #100, Golden, CO, 80403. All results must be received by Wibby Environmental prior to the study closing date shown on the Data Reporting Sheets.

SAFETY

These standards are designed for use by laboratory professionals who are familiar with handling environmental reference materials as well as hazardous materials. If you have any questions about the safe handling of these standards or require a Material Safety Data Sheet (MSDS,) please contact Wibby Environmental at 1-866-WibbyPT (866-942-2978).

QUESTIONS?

If you have any questions regarding these standards or reporting requirements, please call Wibby Environmental at 1-866-WibbyPT (866-942-2978).

6390 Joyce Drive
100
Golden, CO 80403

Phone 303-940-0033
Fax 866-283-0269
www.wibby.com

Soil / Hazardous Waste PT Concentration Ranges and PTRLs

Definitions:

PTRL

NELAC Proficiency Testing Reporting Limits (PTRLs) are provided as guidance to laboratories analyzing NELAC PT samples. At a minimum, the laboratory should use a method that is sensitive enough to generate quantitative results at the PTRLs shown. (REF: NELAC PT FOT Tables)

NA

Not Applicable (NA) has been applied to analytes where a PTRL is not applicable and to state specific analytes that have not had a PTRL determined by the applicable accrediting agency.

TCLP Carbamates (PT-TCLPCARB-SOIL)

NELAC Code	Analyte	Units	Concentration Range	PTRL
Additional State Specific Analytes				
7010	Aldicarb	µg/L	7.50 - 25.0	NA
7015	Aldicarb Sulfone	µg/L	9.50 - 25.0	NA
7020	Aldicarb Sulfoxide	µg/L	7.50 - 25.0	NA
8080	Baygon	µg/L	15.0 - 70.0	NA
7195	Carbaryl	µg/L	10.0 - 50.0	NA
7205	Carbofuran	µg/L	7.50 - 75.0	NA
7710	3-Hydroxycarbofuran	µg/L	7.50 - 37.5	NA
7800	Methiocarb	µg/L	15.0 - 70.0	NA
7805	Methomyl	µg/L	7.50 - 45.0	NA
7940	Oxamyl (Vydate)	µg/L	15.0 - 40.0	NA