

Instructions for: WP Acids (cat # PT-ACIDS-WP)

STANDARD DESCRIPTION

- The WP Acids standard is provided in a flame sealed ampule that contains approximately 2.1 mL of concentrate.
- The WP Acids standard is packaged in acetone.
- Store the standard in the unopened ampule refrigerated (at ~4°C).

ADDITIONAL INFORMATION

- The standard has been provided as a concentrate that must be diluted prior to analysis.
- The standard has been designed and manufactured in compliance with NELAC/EPA criteria. As such each lot of the WP Acids standards will contain a minimum of 80% of the total number of NELAC Accreditation and Experimental Analytes listed on the Data Reporting Sheets and in these instructions.
- Refer to "Reporting Instructions" section of this booklet for guidance on reporting results for analytes that you do not detect.
- The standard should be prepared and analyzed as soon as possible after dilution.

STANDARD PREPARATION, ANALYSIS and STORAGE

1. For best results, the PT standard should be stored refrigerated and then brought to room temperature (near 20°C) when used.
2. Add approximately 990 mL of organic free deionized water to a 1000 mL class A volumetric flask.
3. Carefully open the ampule by snapping off the top at the narrow part of the neck.
4. Transfer exactly 1.0 mL of the PT standard concentrate to the flask using a gas tight syringe and delivering the aliquot below the surface of the water.
5. Bring the volumetric flask to volume with organic free deionized water.
6. Mix the solution by inverting the volumetric flask a minimum of three times.
7. The standard is now ready for preparation and analysis per your routine method(s).
8. The sample should be analyzed as soon as possible after dilution,
9. Report all results in µg/L per the reporting instructions contained in this booklet.
10. Store the diluted standard and any remaining concentrate refrigerated (at ~4°C).



6390 Joyce Drive
#100
Golden, CO 80403

PHONE 303-940-0033
FAX 303-940-0043
www.wibby.com

**Instructions for:
WP Acids (cat # PT-ACIDS-WP) - Continued**

CONCENTRATION RANGE and PTRL

- After preparation per these instructions, the standard will contain the analytes in the following table at a certified Concentration within the ranges shown.
- The NELAC Proficiency Testing Reporting Limit (PTRL) is provided as guidance when analyzing NELAC PT standards. At a minimum, the laboratory should use a method that is sensitive enough to generate quantitative results at the PTRL shown.

Analyte	Units	Concentration Range	PTRL
4-Chloro-3-methylphenol	µg/L	30 - 200	10
2-Chlorophenol	µg/L	30 - 200	10
2,4-Dichlorophenol	µg/L	40 - 190	11
2,4-Dimethylphenol	µg/L	65 - 200	10
2,4-Dinitrophenol	µg/L	100 - 180	10
2-Methyl-4,6-Dinitrophenol	µg/L	60 - 200	14
2-Methylphenol	µg/L	50 - 200	9.5
4-Methylphenol (and/or 3-Methylphenol)	µg/L	50 - 200	5
2-Nitrophenol	µg/L	50 - 190	16
4-Nitrophenol	µg/L	100 - 180	10
Phenol	µg/L	100 - 200	10
Pentachlorophenol	µg/L	55 - 200	11
2,4,5-Trichlorophenol	µg/L	50 - 200	19
2,4,6-Trichlorophenol	µg/L	50 - 200	16
NELAC Experimental Analytes			
2,6-Dichlorophenol	µg/L	40 - 190	11
Additional State Specific Analytes			
Benzoic acid	µg/L	10 - 200	NA
3-Nitrophenol	µg/L	10 - 200	NA
2,3,4,6-Tetrachlorophenol	µg/L	50 - 200	NA