

**Instructions for:  
WP Phthalates and Adipates (cat # PT-PHAD-WP)**

**STANDARD DESCRIPTION**

- The WP Phthalates and Adipates standard is provided in a flame sealed ampule that contains approximately 2.1 mL of concentrate.
- The WP Phthalates and Adipates standard is packaged in acetone.
- Store the standard in the unopened ampule refrigerated (at ~4°C).

**ADDITIONAL INFORMATION**

- The standard has been provided as a concentrate that must be diluted prior to analysis.
- The standard should be prepared and analyzed as soon as possible after dilution.

**STANDARD PREPARATION, ANALYSIS and STORAGE**

1. For best results, the PT standard should be stored refrigerated and then brought to room temperature (near 20°C) when used.
2. Add approximately 990 mL of organic free deionized water to a 1000 mL class A volumetric flask.
3. Carefully open the ampule by snapping off the top at the narrow part of the neck.
4. Transfer exactly 1.0 mL of the PT standard concentrate to the flask using a gas tight syringe and delivering the aliquot below the surface of the water.
5. Bring the volumetric flask to volume with organic free deionized water.
6. Mix the solution by inverting the volumetric flask a minimum of three times.
7. The standard is now ready for preparation and analysis per your routine method(s).s).
8. The sample should be analyzed as soon as possible after dilution,
9. Report all results in µg/L per the reporting instructions contained in this booklet.
10. Store the diluted standard and any remaining concentrate refrigerated (at ~4°C).

6390 Joyce Drive  
# 100  
Golden, CO 80403

Phone 303-940-0033  
Fax 866-283-0269  
www.wibby.com

## Water Pollution Proficiency Testing Concentration Ranges and PTRLs

### Definitions:

#### **PTRL**

NELAC Proficiency Testing Reporting Limits (PTRLs) are provided as guidance to laboratories analyzing NELAC PT samples. At a minimum, the laboratory should use a method that is sensitive enough to generate quantitative results at the PTRLs shown. (REF: NELAC PT FOT Tables)

#### **NA**

Not Applicable (NA) has been applied to analytes where a PTRL is not applicable and to state specific analytes that have not had a PTRL determined by the applicable accrediting agency.

### WP Phthalates and Adipates (PT-PHAD-WP)

NELAC Code	Analyte	Units	Concentration Range	PTRL
<b>Additional State Specific Analytes</b>				
5670	Butyl benzyl phthalate	µg/L	10.0 - 50.0	NA
5925	Di-n-butylphthalate	µg/L	10.0 - 50.0	NA
6062	Di(2-Ethylhexyl) Adipate	µg/L	8.00 - 50.0	NA
6065	Di(2-Ethylhexyl) Phthalate	µg/L	9.00 - 50.0	NA
6070	Diethylphthalate	µg/L	10.0 - 50.0	NA
6135	Dimethylphthalate	µg/L	10.0 - 50.0	NA
6200	Di-n-octylphthalate	µg/L	10.0 - 50.0	NA
6380	1-Methylnaphthalene	µg/L	10.0 - 50.0	NA
6385	2-Methylnaphthalene	µg/L	10.0 - 50.0	NA