

Instructions for: WP Fixed Solids

STANDARD DESCRIPTION

- The WP Fixed Solids standard is provided as a 500 mL ready-to-use sample.
- The WP Fixed Solids standard is unpreserved.
- Store the standard in the unopened bottle at room temperature (~20 - 25°C).

ADDITIONAL INFORMATION

- The WP Fixed Solids standard contains fixed suspended solids. To ensure a uniform sub-sample, shake the bottle vigorously prior to removing an aliquot for the determination of Fixed Suspended Solids.

STANDARD PREPARATION, ANALYSIS and STORAGE

1. For best results, the standard should be stored near room temperature (~ 20 - 25°C).
2. The standard is ready for preparation and analysis per your routine method(s) as received. No dilutions are necessary.
3. Report all results in mg/L per the reporting instructions contained in this booklet.
4. Store any remaining standard at ~ 20 - 25°C.

CONCENTRATION RANGE and PTRL

- After preparation per these instructions, the standard will contain the analytes in the following table at a certified Concentration within the ranges shown.
- The NELAC Proficiency Testing Reporting Limit (PTRL) is provided as guidance when analyzing NELAC PT standards. At a minimum, the laboratory should use a method that is sensitive enough to generate quantitative results at the PTRL shown.

Analyte	Units	Concentration Range	PTRL
Fixed Dissolved Solids	mg/L	140 – 650	NA
Fixed Suspended Solids	mg/L	23.0 – 100	NA
Total Fixed Solids	mg/L	140 - 675	NA